

Product Datasheet

Product code: **AF4535ZTE**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AF4535ZTE offers wireless-ready capability (with the optional WF-20 adapter), a magnification range of 10x - 330x, and a 1.3 Megapixel Edge sensor for detailed imaging. It includes an adjustable polarizer, Automatic Magnification Reading (AMR) for measurement-intensive tasks, and Flexible LED Control (FLC) for tailored lighting options.

Special features: Polarizer - FLC - AMR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-330X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch