

Product Datasheet

Product code: **AF4535ZTLE**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AF4535ZTLE features a 1.3 Megapixel Edge sensor and offers a magnification range of 10x-200x for precise imaging. It is equipped with an adjustable polarizer, Automatic Magnification Reading (AMR) for measurement-intensive applications, and Flexible LED Control (FLC) for customizable lighting. The device is also wireless-ready, requiring the optional WF-20 adapter.

Special features: Polarizer - FLC - AMR - LWD

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-200X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch