

Product Datasheet

Product code: **AF4935ZTE**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AF4935ZTE combines advanced imaging capabilities with features such as Extended Depth of Field (EDOF) and Extended Dynamic Range (EDR), which function exclusively under Windows OS. It offers a magnification range of 10x-330x, a 1.3 Megapixel Edge sensor, and an adjustable polarizer for versatile applications. Additional highlights include Automatic Magnification Reading (AMR), Flexible LED Control (FLC), and wireless readiness with the optional WF-20 adapter.

Special features: Polarizer - FLC - AMR - EDOF - EDR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-330X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch