

Product Datasheet

Product code: **AF7915MZT**

Product type: Digital Microscopes

Made in: Taiwan

Description: Get the detail you need on the go with the Dino-Lite AF7915MZT high-resolution digital microscope. This 5-megapixel handheld digital microscope includes advanced features like EDOF, FLC, AMR, and EDR, offering versatility for accomplishing a wide variety of tasks. AF7915MZT is compatible with the WF-20 Wi-Fi streamer. EDOF, EDR, and DPQ are available on Windows PC only.

Special features: Polarizer - FLC - AMR - EDOF - EDR

Resolution: 5 Megapixel (2592x1944)

Magnification: 10-220X

Picture:



Product dimensions:	OAL	115.0 mm – 4.528 inch
	Width	34.0 mm – 1.339 inch
	Height	34.0 mm – 1.339 inch