

Product Datasheet

Product code: **AF7915MZTL**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AF7915MZTL is part of the Long Working Distance (LWD) range and with its Long Working Distance and magnification range from 10 - 140x, it is suitable for a wide range of applications where extra distance to the object or a larger field of view is required.

Special features: Polarizer - FLC - AMR - EDOF - EDR - LWD

Resolution: 5 Megapixel (2592x1944)

Magnification: 10-140X

Picture:



Product dimensions:	OAL	115.0 mm – 4.528 inch
	Width	34.0 mm – 1.339 inch
	Height	34.0 mm – 1.339 inch