

Product Datasheet

Product code: **AM3111**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AM3111 digital microscope is an economical pick that provides smooth imaging and balanced color performance, making it suitable for a variety of applications, such as simple industrial inspections, educational use, or personal hobbies.

Special features: IR - 60fps - 8 LEDs

Resolution: VGA (640x480)

Magnification: 20-50X, 200X

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch