

DINO-LITE AM4013MT



The Dino-Lite AM4013MT is a convenient, easy to use, professional handheld digital microscope providing superior image quality with up to 200x optical magnification.

It is equipped with a high quality 1.3 Megapixels sensor with a crystal clear image quality and natural color reproduction, enabling you to view and capture objects with the greatest level of detail. This model comes with the DinoCapture 2.0 software for Windows and DinoXcope software for MacOS, allowing you to do calibrated measurements on pictures taken with the Dino-Lite. The aluminum housing offers enhanced protection and durability and has the best look and feel.

Model	AM4013MT
Resolution	1.3 Megapixel (1280x1024)
Magnification	20x - 50x, 200x
Interface	USB 2.0
Nr. of LED's	8
LED on/off switchable	Yes
LED color	White
Housing	Metal housing
Microtouch	Yes
ESD safe	Yes
Output	Picture, video
Video Frame Rate	up to 30 frames-per-second
Measurement	Yes
Calibration	Yes
Diffuser available:	No
Infrared filter:	IR cut-filter >650 nm
Polarizer:	No
Compatibility	Windows 7,8,10 & 11, MacOS 10.9 and up
Software (Included)	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Special feature	No
Package contents	Microscope, Carry pouch, software CD, calibration target, user manual