

Product Datasheet

Product code: **AM4013MT**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite AM4013MT is a convenient, easy to use, professional handheld digital microscope providing superior image quality with up to 200x optical magnification.

Special features: IR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-70X/200X

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch