

## DINO-LITE AM4013MZTL



By using a special polarization filter, the AM4013MZTL is the right choice when working with high glare materials such as plastics and metals. This model features an aluminum alloy housing for enhanced protection and durability and the best look and feel. This model also has a long working distance.

<b>Model</b>	AM4013MZTL
<b>Resolution</b>	1,3 megapixel (1280x1024 pixels)
<b>Magnification</b>	10-90x
<b>Interface</b>	USB 2.0, High speed USB transmission
<b>Nr. of LED's</b>	8
<b>LED on/off switchable</b>	Yes
<b>LED color</b>	White
<b>Housing</b>	Metal housing
<b>Microtouch</b>	Yes
<b>Magnification lock</b>	No
<b>ESD safe</b>	Yes
<b>Output</b>	Picture, video
<b>Video Frame Rate</b>	up to 30 frames-per-second, 15fps at 1.3M
<b>Measurement</b>	Yes (line, radius, circle, 3-point circle, etc.)
<b>Calibration</b>	Yes
<b>Filter/Diffuser</b>	Yes, linear
<b>Compatibility</b>	Windows XP, Vista, 7, 8 & 10, MacOS 10.9 and up
<b>Software (Included)</b>	Windows: DinoCapture 2.0, Mac: DinoXcope
<b>Package contents</b>	Microscope, Carry pouch, Software CD, Calibration target, ZT-Z-CC1 Cap closure, User manual