

Product Datasheet

Product code: **AM4113T**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite AM4113T is a convenient, easy to use, professional handheld digital microscope providing superior image quality with up to 200x optical magnification. It is equipped with a high quality 1.3 Megapixels sensor with a crystal clear image quality and natural color reproduction, enabling you to view and capture objects with the greatest level of detail. This model comes with the DinoCapture 2.0 software for Windows and DinoXcope software for MacOS, allowing you to do calibrated measurements on pictures taken with the Dino-Lite. The AM4113T is the most popular choice for quality inspection in industrial environments, education, office, or in the laboratory.

Special features: IR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-70X, 200X

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch