

Product Datasheet

Product code: **AM4113TL**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AM4113TL is a special model with long working distance up to 15 cm at 20x magnification and around 4cm at 90x magnification. The extra working distance makes the TL series an ideal solution for tasks where you need more working distance to the object, for example repair, rework or assembly tasks.

Special features: IR - LWD

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-90X

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch