

Product Datasheet

Product code: **AM4113ZT**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite AM4113ZT is a convenient, easy to use, professional handheld digital microscope providing superior image quality with up to 200x optical magnification. It is equipped with a high quality 1.3 Megapixels sensor with a crystal clear image quality and natural color reproduction, enabling you to view and capture objects with the greatest level of detail. This model comes with the DinoCapture 2.0 software for Windows and DinoXcope software for MacOS, allowing you to do calibrated measurements on pictures taken with the Dino-Lite. The built-in adjustable polarizer reduces glare and reflection on shiny objects.

Special features: IR - Polarizer

Picture:



Product dimensions:	OAL	115.0 mm – 4.528 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch