

DINO-LITE AM4113ZTL

By using a special polarization filter, the AM4113ZTL is the right choice when working with high glare materials such as plastics and metals. This model also has a long working distance.



Model	AM4113ZTL
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10-90x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	Yes
LED color	White
Housing	Composite/ plastic housing
Microtouch	Yes
Magnification lock	Yes
ESD safe	No
Output	Picture, video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	Yes (line, radius, circle, 3-point circle, etc.)
Calibration	Yes
Filter/Diffuser	Yes, linear
Compatibility	Windows XP, Vista, 7, 8 & 10, MacOS 10.9 and up
Software (Included)	Windows: DinoCapture 2.0, Mac: DinoXcope
Package contents	Microscope, Carry pouch, Software CD, Calibration target, ZT-Z-CC1 Cap closure, User manual