

Product Datasheet

Product code: **AM4115ZTW**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite Edge models provide superior image quality and a high level of flexibility due to the adaptable cap design. The built-in adjustable polarizer reduces glare and reflection on shiny objects. This model can work with 2 magnification rates while maintaining the same working distance, making this the ideal model for macro view applications with a larger field of view and a zoom function to reveal smaller details. Different caps are supplied with the microscope making the Dino-Lite Edge the perfect choice for any professional application.

Special features: IR - Polarizer

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch