

## Product Datasheet

**Product code:** **AM4515ZT**

**Product type:** Digital Microscopes

**Made in:** Taiwan

**Description:** The Dino-Lite Edge models provide superior image quality and a high level of flexibility due to the adaptable cap design. The built-in adjustable polarizer reduces glare and reflection on shiny objects. With the built-in automatic magnification reading, measurements can be performed easily and quickly. Different caps are supplied with the microscope making the Dino-Lite Edge the perfect choice for any professional application.

**Special features:** IR - Polarizer

**Resolution:** 1.3 Megapixel (1280x1024)

**Magnification:** 20-220X

**Picture:**



<b>Product dimensions:</b>	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch