

Product Datasheet

Product code: **AM4515ZT**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite Edge models provide superior image quality and a high level of flexibility due to the adaptable cap design. The built-in adjustable polarizer reduces glare and reflection on shiny objects. With the built-in automatic magnification reading, measurements can be performed easily and quickly. Different caps are supplied with the microscope making the Dino-Lite Edge the perfect choice for any professional application.

Special features: IR - Polarizer

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 20-220X

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	32.0 mm – 1.260 inch
	Height	32.0 mm – 1.260 inch