

Product Datasheet

Product code: **AM73915MZTL**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite AM73915MZTL is a Long Working Distance model that is part of the High Speed range with USB 3.0 connection. It offers superb uncompressed image quality and color reproduction in a robust, compact and appealing metal housing. The Dino-Lite AM73915MZTL is a USB 3.0 model with a magnification range of 10x - 140x and a Long Working Distance of up to 24 cm. This model is member of the DinoLite Edge family with the latest, cutting-edge optics, an unmatched 5 megapixel sensor and several special features such as the Extended Depth of Field (EDOF), Enhanced Dynamic Range (EDR), Automatic Magnification Reading (AMR), Flexible LED Control (FLC), built-in polarization filter, full metal body, extensive measurement functions, calibration, exchangeable front caps and more. *EDOF/EDR only functions under Windows OS

Special features: LWD - IR - Polarizer - FLC - AMR - EDR

Resolution: 5 Megapixel (2592x1944)

Magnification: 10-140X

Picture:



Product dimensions:	OAL	119.0 mm – 4.685 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch