

DINO-LITE AM7515MZTL



With the use of the latest, cutting-edge optics, a brand new 5 megapixel sensor and several special features such as Flexible LED Control (FLC) and full aluminum body, the Dino-Lite AM7515MZTL is a marvel of technology and the best choice for the high-demanding professional. The Dino-Lite AM7515MZTL offers superb image quality and color reproduction in a robust, compact and appealing housing.

The main features of the AM7515MZTL are:

- 5 Megapixel Edge sensor
- Long Working Distance (LWD)
- 10-140x magnification
- Larger Field Of View (FOV)
- Automatic Magnification Reading (AMR)
- Flexible LED Control (FLC)

Model	AM7515MZTL
Resolution	5 megapixel (2592x1944)
Magnification	10-140x
Interface	USB 2.0
Nr. of LED's	8
LED on/off switchable	Yes
LED color	White
Housing	Metal housing
Microtouch	Yes
ESD safe	Yes
Output	Picture, video
Video Frame Rate	up to 30 frames-per-second
Measurement	Yes
Calibration	Yes
Diffuser available:	Yes (N3C-D included)
Infrared filter:	IR cut-filter >650 nm
Polarizer:	Yes, linear
Compatibility	Windows XP, Vista, 7, 8 & 10, MacOS 10.9 and up
Software (Included)	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Special feature	Aut. Magn. reading, Flexible LED Control (FLC)
Package contents	Microscope, Carry pouch, Software CD, Calibration target, User manual, N3C-O- Open cap, N3C-C- Closed cap, N3C-D- Diffuser cap, N3C-Extension cap, N3C-L- Long cap, N3C-S-Side light cap