

## Product Datasheet

<b>Product code:</b>	<b>AM7915MZT</b>
<b>Product type:</b>	Digital Microscopes
<b>Made in:</b>	Taiwan
<b>Description:</b>	With the use of the latest, cutting-edge optics, a brand new 5 megapixel sensor and several special features, the Dino- Lite AM7915MZT is a marvel of technology and the best choice for the high-demanding professional.
<b>Special features:</b>	IR - Polarizer - FLC - AMR - EDOF - EDR
<b>Resolution:</b>	5 Megapixel (2592x1944)
<b>Magnification:</b>	10-220X
<b>Picture:</b>	



<b>Product dimensions:</b>	OAL	115.0 mm – 4.528 inch
	Width	34.0 mm – 1.339 inch
	Height	34.0 mm – 1.339 inch