

DINO-LITE AM8517MT-FUW

Equipped with dual lighting, the AM8517MT-FUW 8-megapixel microscope features 375 nm LEDs for UV applications and white LEDs for high color fidelity imaging. This microscope is suitable for industrial inspections, forensic analysis, counterfeit detection, and other applications where high-quality imaging and UV lighting is required.



Model	AM8517MT-FUW
Resolution	8.0 Megapixel (3840x2160)
Light/ LED type	375 nm UV
Number of LEDs	4 UV + 4 White
Magnification	20x ~ 220x
Interface	USB 2.0
LED on/off switchable	Yes
Microtouch	Yes
ESD safe	No
Output	Picture, video
Video Frame Rate	10 fps at 8.0 MP, MJPEG 25 fps at 2.0 MP, MJPEG
Measurement	Yes
Calibration	Yes
Diffuser available:	No
Infrared filter:	No
Polarizer:	No
Compatibility	Windows 10/11 macOS 10.14 or later
Special feature	Auto and Manual Refocus
Package contents	Microscope, carry pouch, calibration target, user manual, Microscope, carry pouch, calibration target, user manual, N3C-C (Close Cap) N3C-EB (Extended Open Black Cap) N3C-L (Long Cap) N3C-OB (Black Open Cap) N3C-S (Sidelight Cap)