

DINO-LITE AM8517MZT

Capture vibrant and true-to-life images with the 8-megapixel Dino-Lite AM8517MZT digital microscope. Including a polarizer and featuring AMR and eFLC, AM8517MZT provides the versatility to perform a variety of tasks.



Model	AM8517MZT
Resolution	8.0 Megapixel (3840x2160)
Magnification	10x ~ 220x
Interface	USB 2.0
Working distance	142.6 ~ 11.5 mm
Nr. of LED's	8
LED on/off switchable	Yes
LED color	White
Housing	Metal housing
Microtouch	Yes
ESD safe	No
Output	Picture, video
Video Frame Rate	10 fps at 8.0 MP, MJPEG 25 fps at 2.0 MP, MJPEG
Measurement	Yes
Calibration	Yes
Diffuser available:	No
Infrared filter:	No
Polarizer:	No
Compatibility	Windows 10/11
Special feature	Auto and Manual Refocus
Package contents	N3C-C (Close Cap) N3C-D (Diffuser Cap) N3C-D2 (Opal Diffuser Cap) N3C-E (Extended Open Cap) N3C-L (Long Cap) N3C-O (Open Cap) N3C-S (Sidelight Cap)