

DINO-LITE AM8917MZT4

Featuring EDOF for expanding the depth of field, AM8917MZT4 offers 8-megapixels of clarity at 400x ~ 470x and incorporates features like DPQ, refocus mode, AMR, eFLC, and even a polarizer for assisting professionals in high magnification tasks.



Model	AM8917MZT4
Resolution	8.0 Megapixel (3840x2160)
Magnification	400x ~ 470x
Interface	USB 2.0
LED on/off switchable	Yes
LED color	White
Housing	Metal housing
Microtouch	Yes
ESD safe	No
Output	Picture, video
Video Frame Rate	10 fps at 8.0 MP, MJPEG 25 fps at 2.0 MP, MJPEG
Measurement	Yes
Calibration	Yes
Diffuser available:	No
Infrared filter:	No
Polarizer:	No
Compatibility	Windows 10/11 macOS 10.14 or later
Special feature	EDOF, EDR, DPQ, and Auto/Manual Refocus are available on Windows PC only.
Package contents	Microscope, carry pouch, calibration target, user manual, N3C-O- Open cap, N3C-C- Closed cap, N3C-D- diffuser cap, N3C-E- Extension cap, N3C-L- Long cap, N3C-S- Side light cap