

DINO-LITE AM8917MZTL

Capture stunning images with the AM8917MZTL 8-megapixel long working distance digital microscope. Engineered with a comprehensive set of features, including AMR, this model is designed to excel in both imaging and measurement-oriented tasks.



Model	AM8917MZTL
Resolution	8.0 Megapixel (3840x2160)
Magnification	10x ~ 140x
Interface	USB 2.0
Working distance	237.5 ~ 33.9 mm
Nr. of LED's	8
LED on/off switchable	Yes
LED color	White
Housing	Metal housing
Microtouch	Yes
ESD safe	No
Output	Picture, video
Video Frame Rate	10 fps at 8.0 MP, MJPEG 25 fps at 2.0 MP, MJPEG
Measurement	Yes
Calibration	Yes
Diffuser available:	No
Infrared filter:	No
Polarizer:	No
Compatibility	Windows 10/11
Special feature	EDOF, EDR, DPQ, and Auto/Manual Refocus are available on Windows PC only.
Package contents	N3C-C (Close Cap) N3C-D (Diffuser Cap) N3C-D2 (Opal Diffuser Cap) N3C-E (Extended Open Cap) N3C-L (Long Cap) N3C-O (Open Cap) N3C-S (Sidelight Cap))