

Product Datasheet

Product code: **FC-L-MA1**

Product type: Digital Microscopes

Made in: Taiwan

Description: The mirror angle adapter fits onto the models with long working distance: AD4113TL/AD4116TL/AD4112PTL/AD4013MTL/AD4013MZTL/AD7013MTL etc. Mirror angle: 45°

Picture:



Product dimensions:	OAL	68.0 mm – 2.677 inch
	Width	30.0 mm – 1.181 inch
	Height	30.0 mm – 1.181 inch