

## Product Datasheet

**Product code:** **MPTSS2**

**Product type:** Stainless steel spatulas


**Made in:** Pakistan

**Description:** MPTSS2 Stainless Steel spatula is suitable to collect material. With its double tip rounded body, the MPTSS2 is an auxiliary tool for laboratory use. The spoon tips consist in a short, rounded tip and a large, flat, squared tip. Stainless Steel probes are rust-free and equipped with knurled handles which prevent any possible slips.

**Tip description:** short rounded drop-shaped, large flat squared tip

**Body material:** Stainless Steel (S)

**Picture:**



**Product dimensions:**

OAL	180.0 mm – 7.087 inch
Width	14.5 mm – 0.571 inch
Height	5.5 mm – 0.217 inch

**Tip (1) details:**



**Tip (1) dimensions:**

Tip width	13.5 mm – 0.531 inch
Tip thickness	0.5 mm – 0.020 inch

**Tip (2) details:**



**Tip (2) dimensions:**

Tip width	6 mm – 0.236 inch
Tip thickness	1.2 mm – 0.047 inch