

Product Datasheet

Product code: **MPTSS2**

Product type: Stainless steel spatulas


Made in: Pakistan

Description: MPTSS2 Stainless Steel spatula is suitable to collect material. With its double tip rounded body, the MPTSS2 is an auxiliary tool for laboratory use. The spoon tips consist in a short, rounded tip and a large, flat, squared tip. Stainless Steel probes are rust-free and equipped with knurled handles which prevent any possible slips.

Tip description: short rounded drop-shaped, large flat squared tip

Body material: Stainless Steel (S)

Picture:



Product dimensions:

OAL	180.0 mm – 7.087 inch
Width	14.5 mm – 0.571 inch
Height	5.5 mm – 0.217 inch

Tip (1) details:



Tip (1) dimensions:

Tip width	13.5 mm – 0.531 inch
Tip thickness	0.5 mm – 0.020 inch

Tip (2) details:



Tip (2) dimensions:

Tip width	6 mm – 0.236 inch
Tip thickness	1.2 mm – 0.047 inch