

## Product Datasheet

**Product code:** **MPTSS4**

**Product type:** Stainless steel spatulas

**Made in:** Pakistan

**Description:** MPTSS4 Stainless Steel spatula is suitable to collect material. With its double tip rounded body, the MPTSS4 is an auxiliary tool for laboratory use. The spoon tips consist in a long, curved, pointed tip and a long, curved, rounded tip. Stainless Steel probes are rust-free and equipped with knurled handles which prevent any possible slips.

**Tip description:** long curved pointed tip with arrow-shaped, long curved with rounded tip

**Body material:** Stainless Steel (S)

**Picture:**



**Product dimensions:**

OAL	180.0 mm – 7.087 inch
Width	13.0 mm – 0.512 inch
Height	7.5 mm – 0.295 inch

**Tip (1) details:**



**Tip (1) dimensions:**

Tip width	10 mm – 0.394 inch
Tip thickness	0.3 mm – 0.012 inch

**Tip (2) details:**



**Tip (2) dimensions:**

Tip width	0.6 mm – 0.024 inch
Tip thickness	0.2 mm – 0.008 inch