

Product Datasheet

Product code: **NH-M57-NL.9**

Product type: Surgical Scissors

Made in: Germany


Description: The NH-M57-NL are precision Needle holders, ideal for specific work. The line serration provides the handles with extra grip, also preventing any possible slips. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, blunted

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	140.0 mm – 5.512 inch
Width	35.0 mm – 1.378 inch
Height	8.5 mm – 0.335 inch

Tip dimensions:

Tip width	1.50 mm – 0.059 inch
Tip thickness	2.00 mm – 0.079 inch