

Product Datasheet

Product code: **NH-MJ35.9**

Product type: Surgical Scissors

Made in: Germany

Description: The NH-MJ are high precision micro Needle holders, ideal for miniature work. The small size of the tool head is only suitable for holding micro needles. The line serration provides the handles with extra grip, also preventing any possible slips. The locking mechanism on the handles reduces the risk of accidental drops. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, fine, blunted

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	95.0 mm – 3.740 inch
Width	25.0 mm – 0.984 inch
Height	7.0 mm – 0.276 inch

Tip dimensions:

Tip width	0.70 mm – 0.028 inch
Tip thickness	0.90 mm – 0.035 inch