

Product Datasheet

Product code: **NH-SG47-85.B**

Product type: Surgical Scissors

Made in: Pakistan

Description: The NH-SG47 are precision Needle holders, ideal for specific work. The dual-function instrument combines the functions of needle holder and scissors in one instrument. The line serration provides the tips with extra grip, also preventing any possible slips. The locking mechanism on the handles reduces the risk of accidental drops. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, beveled, rhomb serrated

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	120.0 mm – 4.724 inch
Width	65.0 mm – 2.559 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	2.50 mm – 0.098 inch
Tip thickness	1.60 mm – 0.063 inch