

Product Datasheet

Product code: **NH-SG47-TC125.B**

Product type: Surgical Scissors

Made in: Pakistan

Description: The NH-SG47-TC are precision Needle holders, ideal for specific work. The line serration provides the tips with extra grip, also preventing any possible slips. The locking mechanism on the handles reduces the risk of accidental drops. The Tungsten Carbide insert blades are designed to improve holding performance and extend tool life. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, beveled, rhomb serrated

Tip material: Tungsten carbide (TX)

Coating/finishing: Polished

Handle coating: Gold plated (GP)

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	130.0 mm – 5.118 inch
Width	65.0 mm – 2.559 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	2.20 mm – 0.087 inch
Tip thickness	1.20 mm – 0.047 inch