

## Product Datasheet

**Product code:** **RK-06AE**

**Product type:** Digital Microscopes

**Made in:** Taiwan

**Description:** The Dino-Lite RK-06A stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function. Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects.

**Picture:**



<b>Product dimensions:</b>	OAL	220.0 mm – 8.661 inch
	Width	150.0 mm – 5.906 inch
	Height	270.0 mm – 10.630 inch