

Product Datasheet

Product code: **TSD18**

Product type: Solvent Dispenser

Made in: Singapore

Description: TSD18 is an ESD safe safety solvent dispenser with a capacity of 180 millilitres. Ideal for a wide range of applications in the electronic industry as well as for chemistry and medical laboratories. The dispenser prevents evaporation and contamination, while its geometry also prevents unwanted spills and dripping that may endanger the technician.

Tip material: S (AISI 304)

Body material: HDPE

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	56.0 mm – 2.205 inch
	Height	56.0 mm – 2.205 inch