

## Product Datasheet

|                          |   |
|--------------------------|---|
| <b>Product code:</b>     | <b>TSD18AS</b>  |
| <b>Product type:</b>     | Solvent Dispenser   |
| <b>Made in:</b>          | Singapore   |
| <b>Description:</b>      | TSD18AS is an ESD safe safety solvent dispenser with a capacity of 180 millimeters and anti-splash system. Ideal for a wide range of applications in the electronic industry as well as for chemistry and medical laboratories. The dispenser avoids evaporation and contamination, while its geometry prevents spills and dripping that may endanger the technician. |
| <b>Special features:</b> | Anti-splash   |
| <b>Tip material:</b>     | S (AISI 304)  |
| <b>Body material:</b>    | HDPE  |
| <b>Picture:</b>          |   |



|                            |        |                       |
|----------------------------|--------|-----------------------|
| <b>Product dimensions:</b> | OAL    | 105.0 mm – 4.134 inch |
|                            | Width  | 56.0 mm – 2.205 inch  |
|                            | Height | 56.0 mm – 2.205 inch  |