

Product Datasheet

Product code: **WF-20**

Product type: Digital Microscopes

Made in: Taiwan

Description: Coupling with Dino-Lite AF series, the WF-20 is a Wi-Fi streamer adding further versatility, flexibility, and agility to Dino-Lite microscope. The WF-20 allows AF microscope models to be used wirelessly in connection with an iOS (iPhone/iPad) or Android app on any tablet, smartphone or even in the DinoCapture 2.0 software for Windows computers. The WiFi streamer is ideal for field work or presentations and can stream from one Dino-Lite to multiple devices.

Resolution:

Magnification:

Picture:



Product dimensions:	OAL	58.0 mm – 2.283 inch
	Width	36.0 mm – 1.417 inch
	Height	32.0 mm – 1.260 inch