

## Product Datasheet

**Product code:** **WF4115ZT**

**Product type:** Digital Microscopes

**Made in:** Taiwan


**Description:** The Dino-Lite WF4115ZT is a Wireless model that is part of the Universal series. Go from wireless to wired and back in an instant. Work wired via USB or work wireless with the attachable WF-20 Wireless adapter. The WF4115ZT is equipped with premium features such as Flexible LED control (FLC), integrated polarizer, extensive measurement functionality.

**Special features:** WiFi - Polarizer

**Resolution:** 1.3 Megapixel (1280x1024)

**Magnification:** 20-220X

**Picture:**



<b>Product dimensions:</b>	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch