

Product Datasheet

Product code: **WF4115ZTL**

Product type: Digital Microscopes

Made in: Taiwan

Description: The Dino-Lite WF4115ZTL is a Wireless model that is part of the Long Working Distance (LWD) series. Go from wireless to wired and back in an instant. Work wired via USB or work wireless with the attachable WF-20 Wireless adapter. The WF4115ZTL is equipped with premium features such as Flexible LED Control (FLC), integrated polarizer, extensive measurement functionality.

Special features: WiFi - Polarizer - LWD

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-140X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch