

Product Datasheet

Product code: **WF4515ZTL**

Product type: Digital Microscopes

Made in: Taiwan


Description: The Dino-Lite WF4515ZTL is a Wireless model that is part of the Long Working Distance (LWD) series. Go from wireless to wired and back in an instant. Work wired via USB or work wireless with the attachable WF-20 Wireless adapter. The WF4515ZTL is equipped with premium features such as Flexible LED Control (FLC), integrated polarizer, extensive measurement functionality and Automatic Magnification Reading (AMR) for measurement-intensive applications.

Special features: WiFi - Polarizer - LWD - AMR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-140X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch