

Product Datasheet

Product code: **WF4915ZT**

Product type: Digital Microscopes

Made in: Taiwan


Description: The Dino-Lite WF4915ZT is a Wireless model that is part of the Universal series. Go from wireless to wired and back in an instant. Work wired via USB or work wireless with the attachable WF-20 Wireless adapter. The WF4915ZT is equipped with premium features such as Extended Depth of Field (EDOF), Extended Dynamic Range (EDR), Flexible LED Control (FLC), integrated polarizer, extensive measurement functionality and Automatic Magnification Reading (AMR) for measurement-intensive applications.

Special features: WiFi - Polarizer - AMR - EDOF - EDR

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 20-220X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch